

DISPOSABLE I/A HANDPIECES WITH TEXTURED GRIP

Coaxial and bimanual



Soft Tip

Textured
handpiece surface

DISPOSABLE I/A HANDPIECES WITH TEXTURED GRIP

- Textured handpiece surface for comfortable handling / grip
- Compatible with standard luer connector
- Colour coding for each Gauge size



Disposable Coaxial I/A handpieces



23G Coaxial

Compatible with EquiPhaco 1.8mm



21G Coaxial

Compatible with EquiPhaco 2.2mm



21G Coaxial Soft Tip

Compatible with EquiPhaco 2.2mm



20G Coaxial

Compatible with EquiPhaco 2.4mm



Soft Tip I/A

- Soft contact point to fragile tissue and IOL
- Soft tip gives optic feedback to the surgeon
- Soft aspiration port

Disposable Bimanual I/A handpieces



0.8mm



1.1mm



Coaxial

| PRODUCT CODE | SIZE | INCISION SIZE | ASPIRATION CANNULA | ASPIRATION PORT | ASPIRATION HEAD |
|--------------|--------------|---------------|--------------------|-----------------|-------------------------|
| SMS370P | 23G / 0.65mm | 1.8mm | angled 45° | 0.30mm | polished |
| SMS355P | 21G / 0.8mm | 2.2mm | angled 45° | 0.35mm | polished |
| SMS390S | 21G / 0.8mm | 2.2mm | angled 45° | 0.40mm | silicone tip (soft-tip) |
| SMS325P | 20G / 0.9mm | 2.5mm | angled 45° | 0.35mm | polished |

Bimanual

| PRODUCT CODE | SIZE | INCISION SIZE | ASP. CANNULA | ASP. PORT | ASP. HEAD | IRR. CANNULA | IRR. PORT |
|--------------|--------------|---------------|-----------------|-----------|-----------|-----------------|-----------------|
| SBS320A | 23G / 0.65mm | 0.8mm | aeroplane shape | 0.30mm | textured | aeroplane shape | 2 ports, 0.35mm |
| SBS310A | 21G / 0.8mm | 1.1mm | aeroplane shape | 0.35mm | textured | aeroplane shape | 2 ports, 0.45mm |